



Diphenylphosphine

Product number: 312 | Alternative name: Ph₂PH | CAS number: 829-85-6 | C₁₂H₁₁P | Molecular weight: 186.20

Ph₂PH is a clear and colorless air sensitive pyrophoric liquid (b.p. 280°C [536°F]) that should be handled under inert atmosphere at all times. Ph₂PH has been used extensively for the synthesis of phosphine ligands, but it recently has found uses as a base (pK_a = 21.7).¹




DIPHENYLPHOSPHIDE NUCLEOPHILIC ADDITIONS²⁻⁴



R = alkyl, aryl

X = halide, sulfonate

S_N1-ADDITION OF PSEUDO-BENZYLIC AMINES AND ACETATES^{5,6}




X = OAc, NMe₂

71 – 90% yields

R = H, aryl, phosphine


HYDROPHOSPHONYLATION OF ALKENES^{3,7-9}



R = alkyl, aryl, amides, esters, phosphines, silanes, pyridines, sulfides, ethers


61 – 96% yields

SELECTIVE DEMETHYLATION OF ANISOL SUBSTRATES¹⁰




87 – 88% yields

TRANSITION-METAL CATALYZED CARBON-PHOSPHORUS BOND COUPLING¹¹



75% yields

PHOSPHIDE-ACCELERATED DIASTERESELECTIVE CYANATIONS OF KETONES¹²



>99% yield
>96% endo

References

- ¹ Li, J.-N., Liu, L., Fu, Y., and Guo, Q.-X., *Tetrahedron*, **2006**, *62*, 4453.
- ² Kagan, H.B., in *Asymmetric Synthesis*, Morrison, J.D., Editor. **1985**, Academic Press, Inc.: Orlando, FL. p. 1.
- ³ Uriarte, R., Mazanec, T.J., Tau, K.D., and Meek, D.W., *Inorg. Chem.*, **1980**, *19*, 79.
- ⁴ Coote, S.J., Dawson, G.J., Frost, C.G., and Williams, J.M.J., *Synlett*, **1993**, 509.
- ⁵ Hayashi, T., Mise, T., Fukushima, M., Kagotani, M., Nagashima, N., Hamada, Y., Matsumoto, A., Kawakami, S., Monishi, M., Yamamoto, K., and Kumada, M., *Bull. Chem. Soc. Jpn.*, **1980**, *53*, 1138.
- ⁶ Togni, A., Breutel, C., Schnyder, A., Spindler, F., Landert, H., and Tijani, A., *J. Am. Chem. Soc.*, **1994**, *116*, 4062.
- ⁷ King, R.B. and Kapoor, P.N., *J. Am. Chem. Soc.*, **1971**, *93*, 4158.
- ⁸ DuBois, D.L., Myers, W.H., and Meek, D.W., *J. Chem. Soc., Dalton Trans.*, **1975**, 1011.
- ⁹ Shulyupin, M.O., Kazankova, M.A., and Beletskaya, I.P., *Org. Lett.*, **2002**, *4*, 761.
- ¹⁰ Ireland, R.E. and Walba, D.M., *Org. Synth.*, **1977**, *56*, 44.
- ¹¹ Cai, D., Payack, J.F., Bender, D.R., Hughes, D.L., Verhoeven, T.R., and Reider, P.J., *J. Org. Chem.*, **1994**, *59*, 7180.
- ¹² Wilkinson, H.S., Grover, P.T., Vandenbossche, C.P., Bakale, R.P., Bhongle, N.N., Wald, S.A., and Senanayake, C.H., *Org. Lett.*, **2001**, *3*, 553.

FOR MORE INFORMATION

Please call your Regional Customer Service Center today to learn what Entegris can do for you. Visit entegris.com and select the [Contact Us](#) link to find the customer service center nearest you.

TERMS AND CONDITIONS OF SALE

All purchases are subject to Entegris' Terms and Conditions of Sale. To view and print this information, visit entegris.com and select the [Terms & Conditions](#) link in the footer.



Corporate Headquarters

129 Concord Road
Billerica, MA 01821
USA

Customer Service

Tel +1 952 556 4181
Fax +1 952 556 8022
Toll Free 800 394 4083

Entegris®, the Entegris Rings Design®, and other product names are trademarks of Entegris, Inc. as listed on entegris.com/trademarks. All third-party product names, logos, and company names are trademarks or registered trademarks of their respective owners. Use of them does not imply any affiliation, sponsorship, or endorsement by the trademark owner.

©2018-2020 Entegris, Inc. | All rights reserved. | Printed in the USA | 9000-11196TAN-1120