

Large Volume Parenteral and Small Volume Parenteral Filtration

Enables more rigorous process control to improve sterile drug quality

INTRODUCTION

Large volume parenteral (LVP) and small volume parenteral (SVP) are sterile injectable drugs. There are two broadly used methods to produce a sterile drug product: terminal sterilization and aseptic processing of sterilized unit component. For terminal sterilization, in most cases, the product, container, and closure have low bioburden but are not sterile at the time of filling. The product is subjected to a sterilization process in its final container. With aseptic processing, the product is at higher risk of microbial contamination than with terminal sterilization. In an aseptic filling process, the product, containers, and closure are sterilized separately and then brought together under an extremely high-quality environmental condition. It is a well-accepted principle that sterile drugs should be manufactured using aseptic processing only when terminal sterilization is not feasible.

LVPs are sterile drugs packaged in containers labeled as containing more than 100 mL. Most are terminally sterilized, but some can be sterilized by filtration and aseptically processed. SVPs are sterile drugs packaged in containers labeled as containing less than 100 mL. They can be classified into different types including prefilled syringes, ampules, or lyophilized powders. Many SVPs are heat labile requiring sterilizing filtration and aseptic processing.

Large Volume Parenteral




Figure 1. Typical LVP manufacturing process.

TYPICAL LVP/SVP PRODUCTION PROCESS

To minimize risks of microbial, particulate, and pyrogen contamination, sterile product manufacturing requires special processes. Sterilizing filtration is a critical part of any aseptic process, which is removing microorganisms from a fluid stream without adversely affecting product quality. Typical applications include liquid and air/gas sterile filtration.

Multiple factors will affect the effectiveness of any sterilizing filtration process, including the type and number of microorganisms, the properties of the liquid, the filter design and membrane polymer, and the filtration process parameters.

For terminally sterilized products, the bulk solution process should include a filtration step with a microorganism reduction filter or sterilizing grade filter to control bioburden levels and particulates prior to filling final product containers, Figure 1.

If the product cannot be sterilized in the final container, solutions or liquids should be sterilized by filtration through a sterilizing grade filter. Due to the potential additional risks of a sterile filtration process, as compared with other sterilization processes, a second filtration through a sterile sterilizing-grade filter should be taken into consideration, which is also called redundant filtration, Figure 2.

Small Volume Parenteral




Figure 2. Typical SVP manufacturing process.

SEPARATION GOALS

- Prefiltration – removes the colloids and particles contaminants, extends the service life of downstream sterilizing filters.
- Final filtration – sterilizing filtration is a process that can be validated to consistently yield filtrates that are sterile. This is best accomplished by using filters with very low active pharmaceutical ingredient (API) binding, low extractables, and nonpyrogenic. In addition, sterilizing filters must be integrity-tested and either readily sterilizing or supplied presterilized.

LVP/SVP FILTRATION RECOMMENDATIONS

As your trusted partner, Entegris is dedicated to ensuring you have peace of mind throughout every step of your processes. We have filtration expertise and technologies for every filtration step during LVP/SVP production. Our portfolio will help improve sterile drug quality and optimize your production experience.

Application	Recommendation
Prefiltration	Pharmsteri™ PP cartridge series Pharmsteri PES cartridge series
Bioburden reduction	Pharmsteri PES cartridge series
Final sterile filtration	Pharmsteri PES cartridge series
Venting sterile	Pharmsteri PTFE cartridge

For detailed product information contact Entegris.

REGULATORY COMPLIANCE

1. Filters with silicone O-rings or gaskets meet the requirement of USP <88> Plastic Class VI
2. All components meet the FDA Indirect Food Additive requirements cited in 21 CFR 177-182
3. All components meet the (EC) No 1935/2004
4. Manufactured in accordance with an ISO 9001 and ISO 13485 certified quality management systems
5. Meet the particle requirement for SVP set by USP and China Pharmacopeia



Pharmsteri PES/PP cartridge filters.

APPLICATION SUCCESSES

Entegris had helped numerous customers apply filtration to their formulation processes, and enabled them to successfully remove contamination risks and maintain more rigorous process control. To gain peace of mind in your process, contact Entegris.

FOR MORE INFORMATION

Please call your Regional Customer Service Center today to learn what Entegris can do for you. Visit entegris.com and select the [Contact Us](#) link to find the customer service center nearest you.

TERMS AND CONDITIONS OF SALE

All purchases are subject to Entegris' Terms and Conditions of Sale. To view and print this information, visit entegris.com and select the [Terms & Conditions](#) link in the footer.



Corporate Headquarters

129 Concord Road
Billerica, MA 01821
USA

Customer Service

Tel +1 952 556 4181
Fax +1 952 556 8022
Toll Free 800 394 4083

Entegris®, the Entegris Rings Design®, and other product names are trademarks of Entegris, Inc. as listed on entegris.com/trademarks. All third-party product names, logos, and company names are trademarks or registered trademarks of their respective owners. Use of them does not imply any affiliation, sponsorship, or endorsement by the trademark owner.

©2022 Entegris, Inc. | All rights reserved. | Printed in the USA | 4463-12483ENT-0822