AccuSizer® A2000 USP Liquid Particle Counter

The most advanced instrument available for USP <788> Particulate Matter Injections Testing

The new AccuSizer® A2000 USP system is the most advanced instrument available for USP <788> testing. It meets or exceeds all requirements in USP <788> by providing size and count data at the required 10 and 25 μm and easily passes all system standardization tests described in USP <1788>. The unique technology designed into the A2000 USP counter provides capabilities not available in any other liquid particle counter due to each component being the highest specification/performance on the market. The A2000 USP is not just a liquid particle counter; it is a sophisticated particle size analyzer as well.

- Measure size and concentration from $0.5 400 \mu m$
- Sample conservation after measurement
- 64 size channels defining the complete distribution
- 21 CFR 11 features with security management
- Up to 32 user defined channels

Classical light obscuration sensors familiar to pharmaceutical scientists performing USP <788> testing typically have a lower particle size limit near 2 μm . The LE400 sensor uses two detectors (extinction + scattering) to extend the range to 0.5 – 400 μm .

A 64 channel pulse height analyzer provides high resolution results so this is not just a counter – use it as a particle size analyzer on any of your other samples. The AccuSizer software allows conversion from number to volume distribution so results can be compared to other techniques like laser diffraction.

The AccuSizer software is a sophisticated package that automates USP <788> testing including pass/fail criteria and acceptance approval.



AccuSizer A2000 USP System



LE400 sensor principle

Security Settings	х
General	
Max Login Attempts:	3
☐ Enable User Timeout (mins)	15
Password Options	
Minimum Password Length:	6
Password Expiration Interval (days):	30
Special Characters Required	1
☐ Uppercase Required	1
Lowercase Required	1
Numbers Required	1
Ok Cancel	

AccuSizer software security settings




SPECIFICATIONS*


Power

Principle	Single particle optical sizing (SPOS) for high resolution particle size and concentration (particles/mL) analysis. Counts and sizes particles individually, not an ensemble method.
Configurations	Includes sensor, counter, and syringe sampler system. Syringe volumes of 0.5 – 25 mL, 64 size channels, conforms to USP <788>, 21 CFR Part 11 software.
Sensor	LE400-05; $0.5-400~\mu m$ when used alone, light extinction and scattering, summation calibration, particle sensitivity to 10 ppt, concentration limit 10,000 particles/mL, size accurac 2%, count accuracy 10%, recommended flow rate = 30 mL/min, but can be calibrated at othe flow rates depending on configuration.
Sample	$150\;\mu\text{L}-25\;\text{mL}$ (or larger with multiple syringe pulls). Sample is recovered after the measurement process.
Options	Autosampler
	Magnetic stirrer for autosampler
	IQ/OQ documentation for user or complete PSS installation

100 - 120 VAC, 60 Hz or 220 - 240 VAC, 50 Hz

Screen captures from the USP <788> software features





^{*}Sample dependent and may require hardware options, subject to change without notice.

FOR MORE INFORMATION

Please call your Regional Customer Service Center today to learn what Entegris can do for you. Visit <u>entegris.com</u> and select the <u>Contact Us</u> link to find the customer service center nearest you.

TERMS AND CONDITIONS OF SALE

All purchases are subject to Entegris' Terms and Conditions of Sale. To view and print this information, visit entegris.com and select the Terms & Conditions link in the footer.



Corporate Headquarters 129 Concord Road Billerica, MA 01821 **Customer Service**

Tel +1 952 556 4181 Fax +1 952 556 8022 Toll Free 800 394 4083

Entegris®, the Entegris Rings Design®, and other product names are trademarks of Entegris, Inc. as listed on entegris.com/trademarks. All third-party product names, logos, and company names are trademarks or registered trademarks of their respective owners. Use of them does not imply any affiliation, sponsorship, or endorsement by the trademark owner.